

CBC-CAL PLUS

for Genrui
HEMATOLOGY CALIBRATOR

CALIBRATOR

Assay values and ranges provided by Genrui

LOT PLUS1221



2022/01/05

Instrument	Parameter	Assigned Value	Acceptable Range
Genrui KT-6610,KT-6600	WBC x 10 ⁹ /L	9.22	8.92 - 9.52
	RBC x 10 ¹² /L	4.47	4.32 - 4.62
	HGB g/L	134	132 - 136
	MCV fL	96.2	94.2 - 98.2
	PLT x 10 ⁹ /L	234	214 - 254
Genrui KT-6610 VET,VH50	WBC x 10 ⁹ /L	8.99	8.69 - 9.29
	RBC x 10 ¹² /L	4.06	3.91 - 4.21
	HGB g/L	126	124 - 128
	MCV fL	96.3	94.3 - 98.3
	PLT x 10 ⁹ /L	276	256 - 296
Genrui KT-6400,KT-6300 version 2.84 or lower	WBC x 10 ⁹ /L	9.6	9.3 - 9.9
	RBC x 10 ¹² /L	4.22	4.07 - 4.37
	HGB g/L	133	131 - 135
	MCV fL	92.4	90.4 - 94.4
	PLT x 10 ⁹ /L	264	244 - 284
Genrui KT-6400,KT-6300 version 2.85 or higher	WBC x 10 ⁹ /L	9.6	9.3 - 9.9
	RBC x 10 ¹² /L	4.46	4.31 - 4.61
	HGB g/L	135	133 - 137
	MCV fL	93.3	91.3 - 95.3
	PLT x 10 ⁹ /L	262	242 - 282
Genrui New KT-6300,KT-6390,KT-6400	WBC x 10 ⁹ /L	9.91	9.61 - 10.21
	RBC x 10 ¹² /L	4.63	4.48 - 4.78
	HGB g/L	133	131 - 135
	MCV fL	93.7	91.7 - 95.7
	PLT x 10 ⁹ /L	248	228 - 268
Genrui KT-60,KT-40	WBC x 10 ⁹ /L	9.83	9.53 - 10.13
	RBC x 10 ¹² /L	4.62	4.47 - 4.77
	HGB g/L	134	132 - 136
	MCV fL	94.7	92.7 - 96.7
	PLT x 10 ⁹ /L	245	225 - 265
Genrui KT-6300 VET	WBC x 10 ⁹ /L	9.0	8.5 - 9.5
	RBC x 10 ¹² /L	4.57	4.27 - 4.87
	HGB g/L	133	128 - 138
	MCV fL	95.5	92.5 - 98.5
	PLT x 10 ⁹ /L	233	203 - 263

- Notes:**
1. For expiry date please refer to the label on tube. Once opened, please use within 1 week if kept at 2-8°C.
 2. Before using, please gently mix well.
 3. All brands and products are trademarks or registered trademarks of their respective companies.

Genrui Biotech Inc.

Address:4-10F, Building 3, Geya Technology Park, Guangming District,
518106, Shenzhen, P.R.China
Tel: +86 755 26835560 ext: 8026 Fax: +86 755 26678789

CBC-CAL PLUS

for Genrui
HEMATOLOGY CALIBRATOR

CALIBRATOR

Assay values and ranges provided by Genrui

LOT PLUS1221



2022/01/05

Instrument	Parameter	Assigned Value	Acceptable Range
Genrui VH 20,30	WBC x 10 ⁹ /L	8.18	7.68 - 8.68
	RBC x 10 ¹² /L	4.32	4.02 - 4.62
	HGB g/L	119	114 - 124
	MCV fL	95.5	92.5 - 98.5
	PLT x 10 ⁹ /L	266	236 - 296

- Notes:**
1. For expiry date please refer to the label on tube. Once opened, please use within 1 week if kept at 2-8°C.
 2. Before using, please gently mix well.
 3. All brands and products are trademarks or registered trademarks of their respective companies.

Genrui Biotech Inc.

Address: 4-10F, Building 3, Geya Technology Park, Guangming District,
518106, Shenzhen, P.R.China

Tel: +86 755 26835560 ext: 8026 Fax: +86 755 26678789